


Anti-Apob Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04114

Item Number: CM0023089



Introduction

Rabbit polyclonal antibody against human APOB (apoB-100/apoB-48). Validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat APOB. Ideal for lipid metabolism, atherosclerosis, and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	APOB Rabbit pAb
Remarks/Alias	FLDB; FCHL2; LDLQC4; apoB-48; apoB-100
Species	Human
GeneID (Human)	338
GeneID	338
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1890-2190 of human APOB (NP_000375.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04114
Protein Weight	516kDa
Shipping	Ice bag