


Anti-Apob Rabbit Polyclonal Antibody (Wb, If/Icc, Elisa) - P04114

Item Number: CM0018369



Introduction

High-quality Rabbit Polyclonal Antibody against human Apolipoprotein B (ApoB), validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat ApoB with a protein weight of 516 kDa.

[Learn More](#)

Parameter	Value
Product Name	ApoB Rabbit pAb
Remarks/Alias	FLDB, FCHL2, LDLQC4, apoB-48, apoB-100, ApoB
Species	Human
GeneID (Human)	338
GeneID	338
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-330 of human ApoB (NP_000375.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04114
Protein Weight	516 kDa
Shipping	Ice bag