

Anti-Apobec3A Rabbit Polyclonal Antibody For Wb & Elisa - P31941

Item Number: CM0018449



Introduction

Rabbit polyclonal antibody against human APOBEC3A (A3A) for Western blot (WB) and ELISA. Recognizes a 23 kDa DNA-editing enzyme involved in antiviral defense and retrotransposon restriction. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	APOBEC3A Rabbit pAb
Remarks/Alias	A3A; ARP3; PHRBN; bK150C2.1; APOBEC3A
Species	Human
GeneID (Human)	200315
GeneID	200315
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human APOBEC3A (NP_663745.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31941
Protein Weight	23kDa
Shipping	Ice bag