


Anti-Apobec3B Monoclonal Antibody For Wb - Q9Uh17

Item Number: CM0006634



Introduction

High-quality anti-APOBEC3B rabbit monoclonal antibody validated for WB and ELISA. Reacts with human APOBEC3B, a DNA deaminase involved in antiviral innate immunity and retrotransposon restriction. Recombinant immunogen (aa 120-220).

[Learn More](#)

Parameter	Value
Product Name	APOBEC3B Rabbit mAb
Remarks/Alias	A3B; ARP4; ARCD3; PHRBNL; APOBEC1L; bK150C2.2; DJ742C19.2; APOBEC3B
Species	Human
Gene ID (Human)	9582
Gene ID	9582
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UH17
Protein Weight	46kDa
Shipping	Ice bag