


Anti-Apobec3C Polyclonal Antibody For Wb - Q9Nrw3

Item Number: CM0015399



Introduction

Anti-APOBEC3C rabbit polyclonal antibody validated for WB, IF/ICC, IP, and ELISA. Detects human APOBEC3C, a DNA cytidine deaminase that restricts retroviruses and retrotransposons. Suitable for antiviral immunity and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	APOBEC3C Rabbit pAb
Remarks/Alias	A3C; PBI; ARP5; ARDC2; ARDC4; APOBEC1L; bK150C2.3; APOBEC3C
Species	Human
Gene ID (Human)	27350
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human APOBEC3C (NP_055323.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q9NRW3
Protein Weight	23kDa
Shipping	Ice bag