

Anti-Apobec3D Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ak3

Item Number: CM0012837

Introduction

Rabbit polyclonal antibody targeting human APOBEC3D, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying retrovirus restriction and innate immunity. Backed by CamelBio's IVD raw-material supply expertise.

[Learn More](#)

Parameter	Specification
Product Name	APOBEC3D Rabbit pAb
Remarks/Alias	A3D; A3DE; ARP6; APOBEC3E; APOBEC3DE; APOBEC3D
Species	Human
GeneID (Human)	140564
GeneID	140564
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 207-386 of human APOBEC3D (NP_689639.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AK3
Protein Weight	47kDa
Shipping	Ice bag