

Anti-Apobec3G Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Hc16

Item Number: CM0013480



Introduction

Rabbit polyclonal anti-APOBEC3G antibody validated for WB, IHC-P, ELISA. Recognizes the 46kDa cytidine deaminase with antiviral activity; cross-reacts with human, mouse, and rat. Ideal for retroviral restriction research.

[Learn More](#)

Parameter	Value
Product Name	APOBEC3G Rabbit pAb
Remarks/Alias	A3G; ARCD; ARP9; ARP-9; CEM15; CEM-15; MDS019; bK150C2.7; dj494G10.1; APOBEC3G
Species	Human
GeneID (human)	60489
GeneID	60489
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-330 of human APOBEC3G (NP_068594.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC16
Protein Weight	46kDa
Shipping	Ice bag