

Anti-Apobec3G Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hc16

Item Number: CM0017344



Introduction

High-quality rabbit polyclonal anti-APOBEC3G antibody validated in WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat APOBEC3G, a key innate immunity factor inhibiting HIV-1 and other retroviruses.

[Learn More](#)

Parameter	Value
Product Name	APOBEC3G Rabbit pAb
Remarks/Alias	A3G; ARCD; ARP9; ARP-9; CEM15; CEM-15; MDS019; bK150C2.7; dJ494G10.1; APOBEC3G
Species	Human
GeneID (Human)	60489
GeneID	60489
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HC16
Protein Weight	46kDa
Shipping	Ice bag