

## Anti-Apoc1 Rabbit Polyclonal Antibody For Ihc-P And Elisa - P02654

Item Number: CM0027119



### Introduction

Rabbit polyclonal antibody against human Apolipoprotein C-I (APOC1), validated for IHC-P and ELISA applications. Detects human and rat APOC1, a key regulator of plasma lipid metabolism. Suitable for cardiovascular and lipoprotein research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | APOC1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | Apo-CI, ApoC-I, apo-CIB, apoC-IB, APOC1                                                                          |
| Species          | Human                                                                                                            |
| Gene ID          | 341                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of human APOC1 (NP_001636.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | IHC-P, ELISA                                                                                                     |
| Cross Reactivity | Human, Rat                                                                                                       |
| SWISS            | P02654                                                                                                           |
| Protein Weight   | 9kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                          |