

Anti-Apoe Rabbit Polyclonal Antibody For Wb And Elisa - P02649

Item Number: CM0024637



Introduction

Rabbit polyclonal antibody against mouse Apoe, suitable for WB and ELISA. Cross-reacts with human APOE. Ideal for lipid metabolism, cholesterol transport, and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	Apoe Rabbit pAb
Remarks/Alias	Apo-E; Apoe
Species	Mouse
GeneID (human)	348
GeneID	11816
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02649
Protein Weight	36 kDa
Shipping	Ice bag