

Anti-Apoe Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P02649

Item Number: CM0014701



Introduction

Anti-APOE rabbit polyclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat APOE. Ideal for lipid metabolism, Alzheimer's disease, cardiovascular, and neurological research.

[Learn More](#)

Parameter	Value
Product Name	APOE Rabbit pAb
Remarks/Alias	AD2; LPG; APO-E; ApoE4; LDLQ5; APOE
Species	Human
GeneID	348
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 41-307 of human APOE (NP_000032.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02649
Protein Weight	36 kDa
Shipping	Ice bag