

Anti-Apoer2/Lrp8 Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - Q14114

Item Number: CM0015479



Introduction

Rabbit polyclonal antibody targeting human ApoER2/LRP8, validated for WB, IF, IHC, ELISA with cross-reactivity to rat. Ideal for investigating Reelin signaling and viral entry.

[Learn More](#)

Parameter	Value
Product Name	ApoER2/LRP8 Rabbit pAb
Remarks/Alias	MC11; LRP-8; APOER2; HSZ75190; ApoER2/LRP8
Species	Human
GeneID (Human)	7804
GeneID	7804
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900 to the C-terminus of human ApoER2/LRP8 (NP_004622.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q14114
Protein Weight	106kDa
Shipping	Ice bag