


Anti-Apoh Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P02749

Item Number: CM0013546



Introduction

High-quality anti-APOH rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects beta-2-glycoprotein 1 (Apo-H) in human and mouse samples, an essential coagulation regulator.

[Learn More](#)

Parameter	Value
Product Name	APOH Rabbit pAb
Remarks/Alias	BG; B2G1; B2GP1; APOH
Species	Human
GeneID (Human)	350
GeneID	350
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-345 of human APOH (NP_000033.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02749
Protein Weight	38kDa
Shipping	Ice bag