

Anti-Apol1 Polyclonal Antibody For Ihc-P, Elisa - O14791

Item Number: CM0031083



Introduction

Rabbit polyclonal anti-APOL1 antibody validated for IHC-P and ELISA. Reacts with human and rat APOL1 (44 kDa). Essential for studying lipid transport, cholesterol metabolism, and trypanosome lytic factor function.

[Learn More](#)

Parameter	Value
Product Name	APOL1 Rabbit pAb
Remarks/Alias	APOL; APO-L; FSGS4; APOL-I; APOL1
Species	Human
Gene ID (Human)	8542
Gene ID	8542
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 179-398 of human APOL1 (NP_003652.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O14791
Protein Weight	44kDa
Shipping	Ice bag