


Anti-Apolipoprotein A1/Apoa1 Monoclonal Antibody For Wb, If/Icc, Elisa - P02647

Item Number: CM0005464



Introduction

High-quality Anti-Apolipoprotein A1/APOA1 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Ideal for cholesterol metabolism, lipid transport, and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	Apolipoprotein A1/APOA1 Rabbit mAb
Remarks/Alias	HPALP2; apo(a); Apolipoprotein A1/APOA1
Species	Human
GeneID (human)	335
GeneID	335
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-267 of human Apolipoprotein A1/APOA1 (NP_000030.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02647
Protein Weight	28kDa
Shipping	Ice bag