

## Anti-Apolipoprotein A2 Polyclonal Antibody For Wb, If/Icc, Elisa - P02652

Item Number: CM0014641

### Introduction

Rabbit polyclonal antibody targeting human Apolipoprotein A2 (ApoA-II), validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Cross-reacts with mouse and rat. Useful for studying HDL metabolism, lipid transport, and cardiovascular disease.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Apolipoprotein A2 Rabbit pAb                                                                                                         |
| Remarks/Alias    | apoAII; Apo-AII; ApoA-II; Apolipoprotein A2                                                                                          |
| Species          | Human                                                                                                                                |
| GeneID (Human)   | 336                                                                                                                                  |
| GeneID           | 336                                                                                                                                  |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Apolipoprotein A2 (NP_001634.1). |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, IF/ICC, ELISA                                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                    |
| SWISS            | P02652                                                                                                                               |
| Protein Weight   | 17kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                              |