

Anti-Apolipoprotein A2 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P02652

Item Number: CM0009060

Introduction

Recombinant rabbit monoclonal antibody targeting human Apolipoprotein A2 (apoA-II), validated for WB, IF-P, ELISA. Detects human, mouse, rat ApoA2. ApoA2 is a major HDL apolipoprotein involved in lipid transport and metabolism, making this antibody ideal for cardiovascular and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	Apolipoprotein A2 Rabbit mAb
Remarks/Alias	apoAII; Apo-AII; ApoA-II; Apolipoprotein A2
Species	Human
GeneID	336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human Apolipoprotein A2 (NP_001634.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02652
Protein Weight	17 kDa
Shipping	Ice bag