


Anti-Apolipoprotein B (ApoB) Monoclonal Antibody For Wb And Elisa - P04114

Item Number: CM0008118



Introduction

ApoB rabbit monoclonal antibody validated for WB and ELISA with rat cross-reactivity. Recognizes human ApoB (516 kDa), the major LDL protein mediating receptor binding and cholesterol transport. Ideal for cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	ApoB Rabbit mAb
Remarks/Alias	FLDB; FCHL2; LDLQC4; apoB-48; apoB-100; ApoB
Species	Human
Gene ID	338
Immunogen	Synthetic peptide within amino acids 50-150 of human ApoB (P04114)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P04114
Protein Weight	516 kDa
Shipping	Ice bag