


Anti-Apolipoprotein C3 Polyclonal Antibody For Wb, If, Ihc, Elisa - P02656

Item Number: CM0017647



Introduction

Rabbit polyclonal antibody against human Apolipoprotein C3 (APOC3), validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for lipid metabolism and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Apolipoprotein C3 Rabbit pAb
Remarks/Alias	Apo-C3; ApoC-3; APOCIII; Apolipoprotein C3
Species	Human
Gene ID (Human)	345
Gene ID	345
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-99 of human Apolipoprotein C3 (NP_000031.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02656
Protein Weight	11kDa
Shipping	Ice bag