


Anti-Apolipoprotein D (Apod) Polyclonal Antibody For Ihc-P And Elisa - P05090

Item Number: CM0019262



Introduction

Rabbit polyclonal antibody targeting human Apolipoprotein D (APOD), validated for IHC-P and ELISA. Suitable for lipid transport, lipoprotein metabolism, and neurological disease research.

[Learn More](#)

Parameter	Specification
Product Name	APOD Rabbit pAb
Remarks/Alias	APOD
Species	Human
GeneID (Human)	347
GeneID	347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human APOD (NP_001638.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P05090
Protein Weight	21kDa
Shipping	Ice bag