


Anti-Apolipoprotein E (ApoE) Rabbit Polyclonal Antibody For Wb - P08226

Item Number: CM0030477



Introduction

Rabbit polyclonal antibody against mouse Apolipoprotein E (ApoE). Suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	ApoE Rabbit pAb
Remarks/Alias	Apo-E
Species	Mouse
GeneID (Human)	348
GeneID	11816
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of mouse ApoE (NP_033826.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08226
Protein Weight	36 kDa
Shipping	Ice bag