

Anti-Apom Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95445

Item Number: CM0026531



Introduction

ApoM Rabbit pAb validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat ApoM. Recombinant immunogen corresponding to amino acids 22-188 of human ApoM. Suitable for lipid transport research.

[Learn More](#)

Parameter	Details
Product Name	ApoM Rabbit pAb
Remarks/Alias	G3a; NG20; apo-M; HSPC336; ApoM
Species	Human
GeneID (Human)	55937
GeneID	55937
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-188 of human ApoM (NP_061974.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95445
Protein Weight	21kDa
Shipping	Ice bag