

Anti-App (Amyloid-Beta Precursor Protein) Rabbit Monoclonal Antibody [Ko Validated] For Wb, Ip, If, Ihc-P, Elisa - P05067

Item Number: CM0003246



Introduction

KO-validated rabbit monoclonal antibody targeting human APP (aa 671-770). Validated for WB, IP, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for Alzheimer's disease and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] APP Rabbit mAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP; PP
Species	Human
GeneID (Human)	351
GeneID	351
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 671-770 of human APP (P05067).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05067
Protein Weight	87kDa
Shipping	Ice bag