

Anti-App Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05067

Item Number: CM0017078



Introduction

Rabbit polyclonal antibody against human Amyloid-beta precursor protein (APP), validated for WB, IF/ICC, and ELISA. Ideal for neurodegeneration research across human, mouse, and rat samples. SWISS: P05067.

[Learn More](#)

Parameter	Value
Product Name	APP Rabbit pAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP; APP
Species	Human, Mouse, Rat
GeneID (Human)	351
GeneID	351
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human APP (NP_000475.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05067
Protein Weight	87kDa
Shipping	Ice bag