

Anti-App Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - P05067

Item Number: CM0023697



Introduction

High-quality rabbit polyclonal antibody against human APP, KO validated for Western blot and ELISA applications. Recognizes full-length APP and amyloid-beta fragments. Ideal for neuroscience research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] APP Rabbit pAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP; PP
Species	Human
GeneID	351
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Aβ42 (NP_000475.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P05067
Protein Weight	87kDa
Shipping	Ice bag