


Anti-Appbp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q92624

Item Number: CM0022566



Introduction

APPBP2 rabbit polyclonal antibody validated for WB and ELISA in human, mouse, rat. Targets amyloid protein-binding protein 2, a key Cul2 substrate-recognition component in protein ubiquitination. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	APPBP2 Rabbit pAb
Remarks/Alias	PAT1; APP-BP2; HS.84084; APPBP2
Species	Human
GeneID (Human)	10513
GeneID	10513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 386-585 of human APPBP2 (NP_006371.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92624
Protein Weight	67kDa
Shipping	Ice bag