


Anti-APPL1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Ukg1

Item Number: CM0011206



Introduction

APPL1 rabbit monoclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Targets the multifunctional adapter protein APPL1 involved in cell signaling and immune response. Ideal for research on endosomal trafficking and metabolism.

[Learn More](#)

Parameter	Value
Product Name	APPL1 Rabbit mAb
Remarks/Alias	APPL; MODY14; DIP13alpha; APPL1
Species	Human
Gene ID (Human)	26060
Gene ID	26060
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human APPL11 (Q9UKG1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKG1
Protein Weight	80 kDa
Shipping	Ice bag