


Anti-APPL1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukg1

Item Number: CM0030478



Introduction

APPL1 Rabbit Polyclonal Antibody for Western blot and ELISA, cross-reacts with human, mouse, and rat APPL1. Ideal for investigating cell proliferation, metabolism, and insulin signaling. 80 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	APPL1 Rabbit pAb
Remarks/Alias	APPL; MODY14; DIP13alpha; APPL1
Species	Human
GeneID (Human)	26060
GeneID	26060
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human APPL1 (NP_036228.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKG1
Protein Weight	80kDa
Shipping	Ice bag