


Anti-APrt Rabbit Monoclonal Antibody (Kd Validated) For Wb, If/Icc, Elisa - P07741

Item Number: CM0010181



Introduction

Rabbit monoclonal antibody against human APRT, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Targets adenine phosphoribosyltransferase, a key enzyme in purine salvage.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] APRT Rabbit mAb
Remarks/Alias	AMP; APRTD
Species	Human
GeneID (human)	353
GeneID	353
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human APRT (NP_000476.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07741
Protein Weight	20kDa
Shipping	Ice bag