


Anti-Aprt Rabbit Polyclonal Antibody For Wb And Elisa - P07741

Item Number: CM0025922



Introduction

Validated rabbit polyclonal antibody targeting human APRT (adenine phosphoribosyltransferase). Suitable for WB and ELISA; cross-reacts with mouse. Ideal for nucleotide metabolism and disease research.

[Learn More](#)

Parameter	Value
Product Name	APRT Rabbit pAb
Remarks/Alias	AMP; APRTD; APRT
Species	Human
GeneID	353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human APRT (NP_000476.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07741
Protein Weight	20kDa
Shipping	Ice bag