


Anti-APrt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P07741

Item Number: CM0017382



Introduction

Rabbit polyclonal antibody against human APRT, validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat APRT (20 kDa). Ideal for AMP salvage pathway research.

[Learn More](#)

Parameter	Value
Product Name	APRT Rabbit pAb
Remarks/Alias	AMP, APRTD, APRT
Species	Human
Gene ID	353
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human APRT (NP_000476.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07741
Protein Weight	20 kDa
Shipping	Ice bag