

Anti-Aptx Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q7Z2E3

Item Number: CM0029334



Introduction

APTX Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, IP, and ELISA applications. Reacts specifically with human aprataxin (MW 41 kDa). Ideal for DNA repair, base excision repair, and neurodegenerative disease studies.

[Learn More](#)

Parameter	Value
Product Name	APTX Rabbit pAb
Remarks/Alias	AOA; AOA1; AXA1; EAOH; EOAHA; FHA-HIT; APTX
Species	Human
GeneID (Human)	54840
GeneID	54840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 93-342 of human APTX (NP_001182178.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q7Z2E3
Protein Weight	41kDa
Shipping	Ice bag