

Anti-Aqp2 Polyclonal Antibody For Ihc-P, Elisa - P41181

Item Number: CM0015569



Introduction

Rabbit polyclonal antibody targeting human AQP2, validated for IHC-P and ELISA. Cross-reacts with mouse. Detects a 29kDa protein. Ideal for renal water homeostasis and collecting duct research.

[Learn More](#)

Parameter	Value
Product Name	AQP2 Rabbit pAb
Remarks/Alias	NDI2; AQP-CD; WCH-CD; AQP2
Species	Human
Gene ID (Human)	359
Gene ID	359
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 147-271 of human AQP2 (NP_000477.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P41181
Protein Weight	29kDa
Shipping	Ice bag