

Anti-Aqp3 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q92482

Item Number: CM0029804

Introduction

High-quality rabbit polyclonal antibody against human AQP3, validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Targets aquaporin-3, essential for glycerol, water, and H2O2 transport in epithelial and skin research.

[Learn More](#)

Parameter	Value
Product Name	AQP3 Rabbit pAb
Remarks/Alias	GIL; AQP-3; AQP3
Species	Human
GeneID (Human)	360
GeneID	360
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 178-292 of human AQP3 (NP_004916.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92482
Protein Weight	32kDa
Shipping	Ice bag