


Anti-Aqp5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P55064

Item Number: CM0007557



Introduction

Rabbit monoclonal antibody against human AQP5, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for studying water transport, salivary gland function, and skin hydration. MW 28 kDa.

[Learn More](#)

Parameter	Value
Product Name	AQP5 Rabbit mAb
Remarks/Alias	PPKB; AQP-5; AQP5
Species	Human
GeneID (Human)	362
GeneID	362
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 165-265 of human AQP5 (NP_001642.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55064
Protein Weight	28 kDa
Shipping	Ice bag