


Anti-Aqp5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55064

Item Number: CM0025545



Introduction

Rabbit polyclonal antibody to Aquaporin-5 (AQP5), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying water channels in salivary gland function, regulatory volume decrease, and eccrine sweat gland biology.

[Learn More](#)

Parameter	Value
Product Name	AQP5 Rabbit pAb
Remarks/Alias	PPKB; AQP-5; AQP5
Species	Human
Gene ID	362
Immunogen	Synthetic peptide corresponding to amino acids 166-265 of human AQP5 (NP_001642.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55064
Protein Weight	28kDa
Shipping	Ice bag