

Anti-Aqp9 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - O43315

Item Number: CM0025049

Introduction

Rabbit polyclonal antibody against human AQP9, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse. For research on glycerol/water transport, hepatic gluconeogenesis, and arsenic toxicity. - O43315

[Learn More](#)

Parameter	Value
Product Name	AQP9 Rabbit pAb
Remarks/Alias	SSC1, AQP-9, T17287, HsT17287, AQP9
Species	Human
GeneID (Human)	366
GeneID	366
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 211-295 of human AQP9 (NP_066190.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43315
Protein Weight	31kDa
Shipping	Ice bag