

Anti-Aqr Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60306

Item Number: CM0027331



Introduction

Rabbit polyclonal antibody targeting RNA helicase aquarius (AQR/IBP160). Suitable for WB, IHC-P, ELISA in human, mouse, rat. This antibody is ideal for studying pre-mRNA splicing, snoRNP biogenesis, and RNA helicase activity.

[Learn More](#)

Parameter	Value
Product Name	AQR Rabbit pAb
Remarks/Alias	IBP160; fSAP164; AQR
Species	Human
GeneID (Human)	9716
GeneID	9716
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human AQR (NP_055506.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60306
Protein Weight	171kDa
Shipping	Ice bag