

## Anti-Aquaporin-1 (Aqp1) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P29972

Item Number: CM0012288



### Introduction

Recombinant rabbit monoclonal antibody targeting human Aquaporin-1 (AQP1). Validated for WB, IF-P, IHC-P, and ELISA; detects endogenous AQP1 in human, mouse, and rat samples. Ideal for water homeostasis and membrane transport studies.

[Learn More](#)

| Parameter        | Specification                            |
|------------------|------------------------------------------|
| Product Name     | Aquaporin-1 (AQP1) Rabbit mAb            |
| Remarks/Alias    | CO; CHIP28; AQP-CHIP; Aquaporin-1 (AQP1) |
| Species          | Human                                    |
| GeneID (Human)   | 358                                      |
| GeneID           | 358                                      |
| Immunogen        | Synthetic Peptide                        |
| Source           | Rabbit                                   |
| Category         | Monoclonal Antibodies                    |
| Application      | WB, IF-P, IHC-P, ELISA                   |
| Cross Reactivity | Human, Mouse, Rat                        |
| SWISS            | P29972                                   |
| Protein Weight   | 29kDa                                    |
| Shipping         | Ice bag                                  |