

Anti-Aquaporin-1 (Aqp1) Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P29972

Item Number: CM0021322



Introduction

High-quality rabbit polyclonal antibody against human Aquaporin-1 (AQP1), validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for water channel and transport studies.

[Learn More](#)

Parameter	Value
Product Name	Aquaporin-1 (AQP1) Rabbit pAb
Remarks/Alias	CO; CHIP28; AQP-CHIP; Aquaporin-1 (AQP1)
Species	Human
GeneID (Human)	358
GeneID	358
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-269 of human Aquaporin-1 (AQP1) (NP_932766.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29972
Protein Weight	29kDa
Shipping	Ice bag