


Anti-Aquaporin-10 (Aqp10) Rabbit Polyclonal Antibody For Wb And Elisa - Q96Ps8

Item Number: CM0029817



Introduction

High-quality rabbit polyclonal antibody against human Aquaporin-10 (AQP10). Validated for WB and ELISA, reacts with human and mouse. Detects a ~32 kDa band; suitable for glycerol transport and metabolic research.

[Learn More](#)

Parameter	Description
Product Name	AQP10 Rabbit pAb
Remarks/Alias	AQPA_HUMAN; AQP10
Species	Human
GeneID (human)	89872
GeneID	89872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-301 of human AQP10 (NP_536354.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96PS8
Protein Weight	32kDa
Shipping	Ice bag