

Anti-Aquaporin-3 (Aqp3) Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q92482

Item Number: CM0026505

Introduction

Rabbit polyclonal antibody against human Aquaporin-3 (AQP3), validated for IHC-P and ELISA. Detects a 32 kDa glycosylated protein involved in glycerol/water transport; suitable for skin and kidney biomarker research.

[Learn More](#)

Parameter	Value
Product Name	AQP3 Rabbit pAb
Remarks/Alias	GIL; AQP-3; AQP3
Species	Human
GeneID (Human)	360
GeneID	360
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 193-292 of human AQP3 (NP_004916.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q92482
Protein Weight	32kDa
Shipping	Ice bag