

Anti-Aquaporin-4 (Aqp4) Rabbit Monoclonal Antibody For Wb - P55087

Item Number: CM0002138



Introduction

Recombinant anti-Aquaporin-4 (AQP4) rabbit monoclonal antibody for WB, IF-F, IF-P, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for investigating brain water homeostasis, glymphatic transport, and amyloid-beta clearance.

[Learn More](#)

Parameter	Value
Product Name	Aquaporin-4 (AQP4) Rabbit mAb
Remarks/Alias	MIWC; WCH4; hAQP4; Aquaporin-4 (AQP4)
Species	Human
GeneID (human)	361
GeneID	361
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 256-323 of human Aquaporin-4 (AQP4) (NP_001641.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55087
Protein Weight	32KDa/34kDa
Shipping	Ice bag