

Anti-Aquaporin-4 (Aqp4) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P55087

Item Number: CM0027664



Introduction

Rabbit polyclonal antibody against human Aquaporin-4 (AQP4), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Targets a key brain water channel involved in glymphatic clearance and brain water homeostasis.

[Learn More](#)

Parameter	Value
Product Name	Aquaporin-4 (AQP4) Rabbit pAb
Remarks/Alias	MIWC; WCH4; hAQP4; Aquaporin-4 (AQP4)
Species	Human
GeneID (human)	361
GeneID	361
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 256-323 of human Aquaporin-4 (AQP4) (NP_001641.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55087
Protein Weight	35kDa
Shipping	Ice bag