

Anti-Araf Polyclonal Antibody For Wb, If/Icc, Elisa - P10398

Item Number: CM0018735



Introduction

Rabbit polyclonal antibody against human ARAF, validated in WB, IF/ICC, and ELISA. Suitable for detecting ARAF in human, mouse, and rat samples. Ideal for mitogenic signaling and myogenic differentiation research.

[Learn More](#)

Parameter	Value
Product Name	ARAF Rabbit pAb
Remarks/Alias	PKS2; A-RAF; ARAF1; RAFA1; ARAF
Species	Human
GeneID (Human)	369
GeneID	369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 371-606 of human ARAF (NP_001645.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10398
Protein Weight	68kDa
Shipping	Ice bag