

Anti-Arap1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96P48

Item Number: CM0015887



Introduction

Rabbit polyclonal antibody targeting human ARAP1, validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Detects mouse and rat ARAP1. Suitable for actin cytoskeleton and membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	ARAP1 Rabbit pAb
Remarks/Alias	CENTD2; cnt-d2; ARAP1
Species	Human
Gene ID (Human)	116985
Gene ID	116985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 834-1133 of human ARAP1 (NP_001128662.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96P48
Protein Weight	162kDa
Shipping	Ice bag