

Anti-Arc (Harc/Arg3.1) Rabbit Polyclonal Antibody For Wb, Elisa - Q7Lc44

Item Number: CM0030389



Introduction

High-affinity rabbit polyclonal antibody against human ARC/Arg3.1, validated for WB and ELISA. Detects human and mouse ARC; ideal for synaptic plasticity and neuroscience research.

[Learn More](#)

Parameter	Specification
Product Name	ARC Rabbit pAb
Remarks/Alias	hArc; Arg3.1; ARC
Species	Human
GeneID (Human)	23237
GeneID	23237
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ARC (NP_056008.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7LC44
Protein Weight	45kDa
Shipping	Ice bag