

Anti-Arcn1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48444

Item Number: CM0024374

Introduction

Coatomer subunit delta (ARCN1) rabbit polyclonal antibody for WB, IF/ICC, ELISA applications. Validated in human, mouse, rat. Part of the coatomer complex involved in Golgi-to-ER transport.

[Learn More](#)

Parameter	Value
Product Name	ARCN1 Rabbit pAb
Remarks/Alias	COPD; SSMG; SRMMD; ARCN1
Species	Human
GenelD (human)	372
GenelD	372
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human ARCN1 (NP_001646.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48444
Protein Weight	57kDa
Shipping	Ice bag