


Anti-Arf1 Polyclonal Antibody For Wb, Elisa - P84077

Item Number: CM0028731



Introduction

Rabbit polyclonal antibody against human ADP-ribosylation factor 1 (ARF1, P84077). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for studies in Golgi trafficking, synaptic plasticity, and host-pathogen interactions.

[Learn More](#)

Parameter	Value
Product Name	ARF1 Rabbit pAb
Remarks/Alias	PVNH8; ARF1
Species	Human
GeneID (Human)	375
GeneID	375
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-181 of human ARF1 (P84077).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84077
Protein Weight	21kDa
Shipping	Ice bag