


Anti-Arf1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P84077

Item Number: CM0008017



Introduction

ARF1 Rabbit mAb for Western blot, IF/ICC, ELISA. Validated in human, mouse, rat. Targets ADP-ribosylation factor 1, a small GTPase crucial for protein trafficking and Golgi function.

[Learn More](#)

Parameter	Value
Product Name	ARF1 Rabbit mAb
Remarks/Alias	PVNH8; ARF1
Species	Human
Gene ID (Human)	375
Gene ID	375
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 70-150 of human ARF1 (P84077).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84077
Protein Weight	21 kDa
Shipping	Ice bag